

Ex-M13K07

Introduction

Ex-M13K07 is a mutant of M13K07. Two codons for Glu (GAA) were substituted to amber codons (TAG) at the 5' region of gIII of M13K07 helper phage genome by site-directed mutagenesis. Ex-M13K07 has a mutant pIII gene that produces a functional wild-type pIII in suppressing *Escherichia coli* strains but does not make any pIII in non-suppressing *E.coli* strains. Packaging phagemids encoding antibody-pIII fusion in F+ non-suppressing *E.coli* strains with Ex-M13K07 enhanced the display level of antibody fragments on the surfaces of recombinant phage particles resulting in an increase of antigen-binding reactivity >100-fold compared to packaging with M13K07 helper phage. Thus, the Ex-M13K07 provides a system for the efficient enrichment of specific binding antibodies from a phage display library and, thereby, increases the chance of obtaining more diverse antibodies specific for target antigens.

Protocol

Storing the Helper Phage

Recommends storing them at -80 °C with 50% Glycerol. As long as there is a -80 °C stock in the lab, the amplified lab prep can be stored at +4 °C.

Titering the Helper Phage

1. Streak out TG1 onto a 2YT plate. Pick a single colony of TG1, inoculate into 2YT and grow to OD600 = 0.6. Meanwhile prepare a serial dilution of the phage in PBS buffer. Expected number of phage should be about 10¹⁰ pfu/ml; therefore dilute phage to create dilutions 10⁻⁴ to 10⁻⁷.

Note:

1. When preparing serial dilutions of bacteriophage, it is usually much more accurate to dilute 10 µl into 100 to create a 10⁻¹ dilution or 10 µl into 1000 µl to create a 10⁻² dilution rather than to dilute 1 into 10 µl or 1 µl into 100 µl because one tends to obtain grossly exaggerated titers resulting from small pipetting errors. This is particularly critical for the first dilution.

2. Add 200 µl of TG1 cells at OD600 = 0.6 to individual sterile culture tubes (the number of phage dilutions you plan to plate) in a test tube rack.

3. Add 100 µl of each serial dilution of helper phage to the culture tubes containing the TG1 cells.

4. Place the test tube rack into a 37 °C water bath for 15 minutes to allow the helper phage to attach to the cells.

5. Meanwhile melt NZY Top Agarose in the microwave and allow it to cool to ~48/50 °C.

You need to pay close attention to the bottle containing the top agarose because it can quickly boil over; alternatively use a low power setting on the microwave.

6. At the conclusion of the 15-minute incubation, add ~3 ml of the NZY top agarose to the test tube containing the helper phage & *E.coli*. Remove the tube from the rack (which is still sitting in the water bath) and give it a quick

flick of your wrist mixing the contents and immediately pour it onto an NZY plate. (I usually do three plates at one by pipetting 10 ml of the NZY top agarose from the bottle, dispensing ~3.3 to one tube, ~3.3 to second tube and then another ~3.3 to a third tube. I then pour the plates starting with the first tube I put the NZY into.)

7. Allow the top agarose to cool (~5 minutes).

8. Invert the plate and incubate overnight at 37 °C.

9. Count the number of plaques.

Amplifying the Helper Phage

1. Streak out TG1 onto a 2YT plate.

2. Pick a single colony of TG1 and inoculate into 5ml 2×YT and grow until OD600 = 0.3 (~2.5×10⁸ cells/ml).

3. Add helper phage at a multiplicity of infection (MOI) of 20:1 (phage-to-cells ratio).

4. Grow the culture at 37 °C with vigorous aeration (~300rpm) for more 8 hours.

Note: (Add kanamycin to a final concentration of 50 µg/ml to the medium 6h after the helper phage and cells have been allowed to grow together.)

5. Add every 1mL of the culture to 100mL 2YT, and culture at 30 °C with 230rpm overnight.

6. Spin culture at 10,000rpm for 20 minutes. Transfer supernatant to a new tube.

7. Add 1/5 volume of 2.5 M NaCl/20% PEG-8000 solution to tube, and let sit at 4 °C for 2h after mixing.

8. Spin at 10,000rpm for 20 minutes. Decant the supernatant, and spin with 1mL 1×PBS briefly.

9. Add 1/5 volume of 2.5 M NaCl/20% PEG-8000 solution to tube, and let sit at 4 °C for 1h after mixing.

10. Spin at 10,000rpm for 20 minutes. Decant the supernatant, and spin with 200 µL 1×PBS again briefly.

11. Add 100% glycerin to 50% concentration and store in the refrigerator at -80 °C.

Related products

Product name	Size	Cat#
pComb3Xss	10 µg	P001
pComb3XTT	10 µg	P002
pComb3XLambda	10 µg	P003
pComb8	10 µg	P004
pCANTAB5E	10 µg	P005
Helper phage M13K07	1 mL	P006
Helper phage VCSM13	1 mL	P007
<i>E.coli</i> TG1	1 mL	P008
<i>E.coli</i> ER2738	1 mL	P009
<i>E.coli</i> XL1-blue	1 mL	P010
<i>E.coli</i> SS320	1 mL	P012
<i>E.coli</i> TOP10 F'	1 mL	P016

This product is for research use only and is not approved for use in humans or in clinical