



Anti-Human IgG, AlpSdAbs[®] VHH(pH-Green 500)

Summary

Code	023-101-011
Immunogen	Recombinant Fc region of human IgG
Host	Alpaca pacous
Isotype	VHH domain of alpaca IgG2b/2c
Conjugate	pH-sensitive Green 500(Ex: 445nm, Em: 503nm), 2 moles Green 500 per mole VHH
Specificity	Human IgG(Fc _γ fragment specific)
Cross-Reactivity	Recognizes human IgG Fc _γ fragment specifically, and reacts with cynomolgus IgG. No Cross-reactivity to rabbit , mouse, rat, goat IgG
Purity	Recombinant Expression and Affinity purified
Concentration	1mg/ml
Formation	Liquid, 10mM PBS (pH 7.4), 50% Glycerol
Storage	Store at -20 °C(Avoid freeze / thaw cycles), Protect from light

Description

Anti-Human IgG, AlpSdAbs[®] VHH(pH-Green 500) is designed for studying on the internalization of antibodies. Anti-Human IgG, AlpSdAbs[®] VHH(pH-Green 500) is based on recombinant single domain antibody to human IgG Fc coupled to pH-Green 500. Based on immunoelectrophoresis and/or ELISA, Anti-Human IgG, AlpSdAbs[®] VHH(pH-Green 500) reacts with the Fc fragment of human IgG heavy chain but not with the Fab portion of human IgG. Anti-Human IgG, AlpSdAbs[®] VHH(pH-Green 500) is an effective detection tool and can be used as a useful tool for the evaluation of antibody potency prior to ADCs.

Background

pH-Green 500 has pH-sensitive fluorescence excitation/emission spectra of 445/503 nm, and its fluorescence emission increases in intensity with increasing acidity. This increase is particularly dramatic in the range pH 4.5–9, as commonly seen within endocytic vesicles. pH-Green 500 is essentially dark in the extracellular environment; however, upon internalization, it elicits a bright fluorescent signal in the acidic environment of the endosomes.

VHH are single-domain antibodies derived from the variable regions of heavy chain of Camelidae immunoglobulin. The size of VHH is extremely small(<15KDa) compared to other forms of antibody fragment, which significantly increase the permeability of VHH. Thus VHH is considered of great value for research, diagnostics and therapeutics.

Benefits

High lot-to-lot consistency
Increased sensitivity and higher affinity
Animal-free production

Application notes

Antibody Internalization Test: 2ug per 10ug test antibody

Dilution factors are presented in the form of a range because the optimal dilution is a function of many factors, such as antigen density, permeability, etc. The actual dilution used must be determined empirically.

This product is for research use only and is not approved for use in humans or in clinical