



Anti-HA tag, AlpHcAbs[®] Mouse IgG2a antibody (HRP)

Summary

Code 003-301-005

Immunogen HA tag fused KLH

Host Alpaca pacous

Isotype VHH domain of alpaca IgG2b/2c fused to Mouse IgG2a Fc(mutation)

Conjugate HRP

Specificity HA tag(YPYDVPDYA)

Cross-Reactivity Highly selective for HA tag sequence

Purity Recombinant Expression and Affinity purified

Concentration 1mg/ml

Formation Liquid, 10mM PBS (pH 7.5), 0.05% sucrose, 0.1% trehalose, 0.01% proclin300,50% Glycerol

Storage Store at -20 °C, (Avoid freeze / thaw cycles), protect from light

Description

Anti-HA tag, AlpHcAbs® Mouse IgG2a antibody(HRP) is designed for detecting HA tag fusion protein specifically. Anti-HA tag, AlpHcAbs® Mouse IgG2a antibody(HRP) is based on recombinant, mouse IgG2a Fc fused single domain antibody to HA tag coupled to HRP. Based on western blot and ELISA, Anti-HA tag, AlpHcAbs® Mouse IgG2a antibody(HRP) reacts with the HA tag selectively, no reactivity with other proteins.

Background

The HA tag is widely used for detecting, manipulating or purifying proteins. This peptide can be expressed and detected with the protein of interest as an amino-terminal or carboxy-terminal fusion. Because of its small size, HA tag is unlikely to affect the tagged protein's biochemical properties. HA tag is useful for the labeling and detection of proteins using immunoblotting, immunoprecipitation, and immunostaining techniques.

Using antibody with Fc(mutation), the background from Fc receptors will be eliminated.

Benefits

High lot-to-lot consistency Increased sensitivity and higher affinity Animal-free production

Suggested Working Concentration

ELISA 1:5000 -1:20000 WB 1:5000 -1:20000

Dilution factors are presented in the form of a range because the optimal dilution is a function of many factors, such as antigen density, permeability, etc. The actual dilution used must be determined empirically.

This product is for research use only and is not approved for use in humans or in clinical

Website: alpvhhs.com E-mail: service@nb-biolab.com Phone: 400-166-9953